
## Dr. Sabharwal's Wound Care




#### Lidocaine 4% + Menthol 1% Patch


#### Licaine +

For temporarily relief of minor pain associated with: operation, arthritis, simple backache, bursitis, tendonitis, muscle strains, sprains & bruises. Desensitize aggravated nerves. Apply for up to 8 hours. Easy application. Individually packed.

#### Lidocaine 4% Patch Licaine 4%

For temporarily relief of minor pain associated with: operation, arthritis, simple backache, bursitis, tendonitis, muscle strains, sprains & bruises. Desensitize aggravated nerves. Apply for up to 8 hours. Easy application. Individually packed.





#### Lidocaine 5% Gel Patch Licaine 350 / 700

Contains: Lidocaine 700mg / 1g Adhesive Lidocaine 350mg / 1g Adhesive

#### Pain Relieving Patch Foot Patch

It is fast acting relief for the most sensitive part of the body containing more the most nerves. The only patch designed with your feet in mind. Camphor, Menthol & Methyl Salicylate deeply penetrate to numb away pain and is maximum amount of pain relief available with out a prescription.




# anarwal's Wounder

#### Dr. Sabharwal's Wound Care

**Factory:** 

137, Buranwala, Barotiwala-174103, HP, India

Telephone: +91 9056770088
Telefax: +91 172 2640019
E-mail: buyers@drsabharwal.in
Website: www.drsabharwal.com



#### Capsicum Plaster Panther Plast C

It contains Capsicum Oleoresin and Methyl Salicylate. Capsicum Plaster is used as a counter-irritant and to produce heat on undamaged skin. It is used to relieve pain in Rheumatic conditions, Remove the patch if heat is unbearable.

#### Ketoprofen 30mg Patch Ketofen

It is effective in a variety of acute and chronic muscular pains.




#### Diclofenac 200mg Patch Panther Plast 200

Diclofenac transdermal patch used to relieve pain caused by sprains, strains, or bruises. Diclofenac is a strong analgesic.



Contains Capsicum Oleoresin 0.5% w/w, Methyl Salicylate 4% w/w, Menthol 0.5% w/w, Camphor 1% w/w, Thymol 1% w/w & Eucalyptus Oil 3% w/w.

A powerful rufificient good for muscular pains or strain, stiff neck & aching shoulders, back ache or lumbago.




#### Belladonna Plaster XL Bella Plast

Belladonna Plaster consists of a fabric coated with adhesive mass containing Belladonna Extract. The plaster is porous & the adhesive surface is protected with printed /plain silicon paper. Provides temporary relief from pain by worming where applied. It is used as counter irritant on undamaged skin. It is good for boil.



Capsicum Plaster consists of a suitable fabric evenly spread with self adhesive mass containing Capsicum Oleoresin. The plaster is porous & the adhesive surface is protected with silicon based release paper.

Capsicum Plaster is used as a counter-irritant and to produce heat on undamaged skin. It is used in the treatment of Rheumatic conditions.




## Motion Sickness Patch Trawell Patch

Patch Contains Promethazine Hydrochloride I.P. 5% w/w in adhesive base.

It is used to prevent or treat motion sickness. It is also used to prevent or treat nausea and vomiting.







#### Adhesive First Aid Dressing of Different Sizes, Shapes & Colour


**Panther Plast FAD** 


Breathable Film, Superior Strong Hypoallergenic Adhesive. Latex Free. Medicated with superior antibacterial Benzalkonium Chloride 0.5% w/w

First Aid remedy for small cuts, wounds & abrasions.

#### Sizes Available:-

38 x 50mm 38 x 38mm 60 x 100mm Oblona Butter Fly Finger Tip 60 x 110mm 19 x 63mm 50 x 76mm Round - 25mm 38 x 60mm 100 x 150mm 19 x 72mm Square 50 x 76mm etc. 100 x 200mm 25 x 76mm 100 x 250mm










#### 20 Multi Colour Strips of 1.9cm x 7.2cm

## Cannula Fixator CanuFix

Used to cover and protect catheter sites to secure devices to the skin. Available in different shape and materials like Transparent Film, Non-Wowoven Substrate and Elastic Adhesive Bandage.









### Burn Gel Dressing Excel Gel

Leno weave cotton gauze impregnated with an ointment of soft paraffin having Chlorhexidine Acetate B.P. 0.5% w/w, Gamma Sterile.



Leno gauze impregnated with an ointment of soft paraffin containing Chlorhexidine Acetate B.P. 0.5% w/w, Gamma Sterile.







Povidone Iodine 10% Tulle Pipdine 10% Tulle

It consists of a knitted synthetic close weave fabric impregnated with water soluble polyethylene glycol (PEG) base containing 10% Povidonelodine.

Neomycin, Bacitracin & Polymixin B Sulphate Tulle
N.T. Biotic Tulle

Tulle Gauze impregnated with ointment containing per 1 gram: Polymixin B Sulphate USP 5000 Units, Bacitracin Zinc USP 400 Units & Neomycin Sulphate USP 3400 Units





## Paraffin Gauze Dressing BP'93 Excel Tulle

Leno gauze impregnated with an ointment of soft paraffin containing, Gamma Sterile.

## Silver Sulphadiazine Gauze Dressing Silva Tulle

Tulle Gauze impregnated with ointment containing Silver Sulphadiazine 1% w/w





#### Poly Urethane Casting Tape/Bandage

#### **Excel Cast**

Polyester based fabric coated with water curable polyurethane. It is of variable setting time. It is light weight & porous. It is X-ray transparent. Patient can return to work early. It is more pliable & soft made especially for tropical climate. It is of high strength & is more conformable. To keep it fresh it has been packed in two shells of alluminium foil & one card board carton.

### Adhesive Tape USP HealoPlast

Cloth coated with hypoallergenic latex free adhesive mass which does not offset when unrolled.



## Plaster of Paris Bandage BP **Healo Cast**

Leno weave fabric coated with plaster of Paris. Low powder loss when dipped in water. The more sticky powder adheres the layers more with each other & makes all layers in one & this one layered cast has more than double the strength as compared to our competitors. Therefore less bandages are required to make a plaster & the resultant plaster cast in of light weight & high strength.

## Cohesive Bandage **Excel Co Bandage**

It is used to keep the dressing at its place or to splint two fingers or toes together. Available in skin colour or 7 colours of rainbow.



plaster of Pa

Bandage B.P.

Plaster of Paris

Bandage B.P.88

### Bandage Skin Traction **Excel Traction**

Mainly composed of adhesive tape having nonallergenic adhesive covered with release paper plastic plate, foam, cord and cotton stretch bandage for giving traction in fracture femur or slipped disk etc. by pulling legs.







## Cotton Crepe Bandage BP Excel Crepe

It is woven on imported computerized knitting machines. The bandage can be of pure cotton, pure nylon or any other material as per customers demand. It is available in small, medium or large size. Patient can choose the best befitting size for him.

#### Elastic Adhesive Bandage BP Excel Plast


Good aesthetic appeal. Thick fabric, allows free muscle & joint movement, Contains natural rubber & wool fat, Porous adhesive to allow skin breathing. Adhesive covered with release liner for easy unrolling,





#### Tubular Net Bandage Excel Net

This is a tubular net bandage which is highly elastic both in length & width wise. Net bandage is indicated to hold the dressing/medication in its place.



# Bandage Elasticated Elasticated cost in ribbed bandage with reser threads

## Elasticated Tubular Bandage Excel Grip

Tubular support bandage is easy and quick to wear, provides even pressure, is comfortable & allows freedom of movement. It is washable & stays in place without pins or tape. It is used for the treatment & prevention of varicose veins, leg oedema, in minor musculine injuries or for joint support. Available in 5cm to 30cm dia.



### Cotton Stockinette **Excel Sleeve**

It is a tubular bandage. It is worn on the limb to give cushioning before plastering. It absorbs sweat & protect bony prominence from plaster cast.

## Knee Support Bandage Kozy Knee Cap

It is woven on imported computerized knitting machines. The bandage can be of pure cotton, pure nylon or any other material as per customers demand. It is available in small, medium or large size. Patient can choose the best befitting size for him.





## Alcohol Wipe **Excel Wipe**

Non-woven 2.5x5cm pad saturated with Isopropyl Alcohol 70% v/v.

For skin disinfection prior to injection/incision.

## Chlorhexide & Isopropyl Wipe Hexi Wipe-C


Wipe saturated with Chlorhexidine Gluconate BP 2% & Isopropyl Alcohol 70%





### Microporous Surgical Tape Excel Pore +

A latex-free, hypoallergenic non-woven paper tape that is gentle to the skin yet adheres well and leaves minimal adhesive residue upon removal. Breathable surgical tape. Available in white or skin colour.



#### Hypoallergenic Surgical Tape **Excel Hypoplast**

Porous Fabric coated with Breathable, Latex Free & Hypoallergenic Adhesive.

#### Corn Removal Pad Corn-Rid

Contains Salicylic Acid 20% to provide relief and remove corns with extra care and maximum comfort.

# Dr Sabharwal's Extension Plaster WETHANGE SENTENDED COLOR SENTENDED CO

## Taffeta Self Adhesive Plaster **Excel Silk Tape**

It consists of a Silk cloth evenly spread with latex free, hypoallergenic self-adhesive mass, It is available in various colours,

#### **Extension Plaster**

Width-wise stretchable fabric with adhesive to provide support to joint without hampering mobility.

Very useful for applying skin traction in fracture femur. Being porous does not cause maceration of skin. Available with Hypoallergenic adhesive too,







#### **Excel Eye Pad**

Highly absorbent cotton wrapped in absorbent nonwoven fabric and cut into oval shape, packed individually, supplied both Sterile & Non-Sterile,

Eye Pad with Adhesive Perfect Eye Pad

An oval shape with Non-Adherent Pad in the centre and with adhesive border, packed individually, supplied both in Sterile & Non-Sterile.





## Nonwoven Surgical Pad **Excel Pad NW**

Highly absorbent cotton wrapped in absorbent non-woven fabric, packed individually, Gamma Sterile.



## Non-Adherent Surgical Pad Sticknil Pad

Thick absorbent pad with non adherent perforated film/net on one or both sides, Gamma Sterile.

# GAUZE SWAB B.P.'888 XL SWAB Store flat in a well ventilated place away from sunlight. Do not accept or use if package is damaged or opened. Sterifized by Gamma Radiation. Do not resterifize the content. FOR EXTERNAL SINGLE USE ONLY. Size: 10cm x 10cm x 12PLY Mid. in India by: D.L. No. MNB/06/348 Dr. Sabharwal's Wound Care 137, Buranwale, Barctiwale, Diett. Solan (H-P)- 174 103 Oustomer Care Tele Fax No. - 91 172 2565789 EGREP J. Sun Jr. 28, 70ch park food, flusation, 80 94 NT UX.

### Gauze Swab B.P. Excel Swab

It is made of absorbent cotton gauze of B.P. Quality type 13 light with all cut threads folded inside. It is made on fully automatic machines. It is used to clean or cover the operated wound or any other wounds after applying antiseptic ointment.

